

## BIOBASE BKM-Z80BIII B Series Table Top Autoclave, New



Capacity 2.8 Cubic Feet, One year warranty.

Rating: Not Rated Yet

**Price**

Sales price \$12,495.00

Tax amount

[Ask a question about this product](#)

Manufacturer [BIOBASE](#)

### Description

- LCD display.
- Internal water filter.
- Water quality monitor.
- Door interlock for safety.
- Design temperature 144° C.
- Temperature range 105° C > 136° C.
- Vacuum system operates quietly.
- Capacity 2.8 Cubic Feet (80 Liters).
- Build-in easy to clean open type water tank.
- Medical-grade 304 Stainless Steel Chamber.
- Dual chamber overpressure protection for safety.
- 208 VAC, 3-Phase power required.
- Date of Manufacture: October 31, 2017.

