


BIOBASE BKM-Z80BIII B Series Table Top Autoclave, New



Capacity 2.8 Cubic Feet, One year warranty. BKMB(III) series sterilizer is an automatic high temperature and pressure rapid sterilizer which works with steam as medium.

Rating: Not Rated Yet

Price

Sales price \$12,495.00

Tax amount

[Ask a question about this product](#)

Manufacturer [BIOBASE](#)

Description

- LCD display.
- Internal water filter.
- Water quality monitor.
- Door interlock for safety.
- Design temperature 144° C.
- Temperature range 105° C > 136° C.
- Vacuum system operates quietly.
- Capacity 2.8 Cubic Feet (80 Liters).
- Build-in easy to clean open type water tank.
- Medical-grade 304 Stainless Steel Chamber.
- Dual chamber overpressure protection for safety.
- 208 VAC, 3-Phase power required.
- Date of Manufacture: October 31, 2017.

