

## Forma Scientific Model 3050 Gas Guard Monitor\Controller



For CO<sub>2</sub> and N<sub>2</sub> Gas

Rating: Not Rated Yet

**Price**

Sales price \$115.00

Tax amount

[Ask a question about this product](#)

Manufacturer [Forma Scientific](#)

### Description

Forma Scientific Model 3050 Gas Guard Monitor\Controller:

- For CO<sub>2</sub> and N<sub>2</sub> Gases
- Hooks up to two (2) gas tanks
- Includes 90-day warranty
- Requires 115 VAC, 50/60 Hz, 0.3 A (single phase) main power
- When one tank goes empty, the Gas Guard goes into audible alarm and switches to secondary tank