

HP 6890 Series GC System



Highly effective for detection of volatile compounds

Rating: Not Rated Yet

Price

Sales price \$37,000.00

Tax amount

[Ask a question about this product](#)

Manufacturer [HP](#)

Description

Agilent 6890 GC / 5973 MS Turbo

Includes dual split/splitless injector, TCD, FID & Mass Spectrometer detectors

NIST Library loaded on HP Compac computer running Chemstation. System given PM before shipping

CTC Analytics Combi Pal Headspace/ Liquid Autosampler includes:

(Allows user to use both liquid and headspace injections)

- Liquid syringe kit
- Agitator/ incubator
- Tray set vial holders to be determined
- Wash station
- Software integrating with Chemstation