

NuAire NU-407FM-600 6 ' Biological Safety Cabinet with Stand



Class II Type A2

Rating: Not Rated Yet

Price

Sales price \$2,695.00

Tax amount

[Ask a question about this product](#)

Manufacturer [Nuair](#)

Description

NuAire NU-407FM-600 6' Biological Safety Cabinet with Stand:

- Provides Class II Type A2 laminar airflow
- Includes UV Light
- Work access opening: 8"
- Internal Work Area: 70.5" W x 24.5" D x 28.0" H
- External Dimensions: 75.0" W x 33.0" D x 83.0" H
- Operates on 115 VAC, 60 Hz main power

Units in box: 1