

## NuAire NU-425-600 6' Biological Safety Cabinet with Stand



Class II Type A2

Rating: Not Rated Yet

**Price**

Sales price \$2,895.00

Tax amount

[Ask a question about this product](#)

Manufacturer [Nuair](#)

### Description

- Class II Type A2.
- Stand is included.
- Manufactured: 2001.
- Ultra quiet while in operation.
- DownFlow Velocity: 70 FPM.
- Inflow Velocity: 105 FPM (.53 m/s).
- Front sash is hinged (10" sash opening).
- Both Fluorescent and UV lights included.
- 110-115 VAC, 60 Hz, 12 Amp circuit required.
- Exterior Dimensions 53.0" W x 33.0" D x 60.0" H.
- Interior Dimensions: 46.0" W x 23.5" D x 28.5" H.