

## Sakura Tissue-Tek Glas Automated Coverslipper



400 slides per hour

Rating: Not Rated Yet

**Price**

Sales price \$2,995.00

Tax amount

[Ask a question about this product](#)

Manufacturer [Sakura](#)

### Description

Sakura Tissue-Tek Glas 6400 Automated Coverslipper

- Benchtop configuration
- Application: Histology, cytology specimens
- Able to prepare up to 400 cover slides per hour
- Load up to 60 slides at a time with three slide trays
- Includes two slide trays
- Volume of the mounting medium is adjustable
- Able to change the coverglass size from 40, 50, 55 or 60 mm
- External Dimensions: 20.0" W x 25.0" D x 25.1" L (50.6 x 63.4 x 63.7 cm)
- Weight: 152 lbs (69 kg)
- Operates on 115 VAC, 50/60 Hz
- Like new condition, well maintained

Units in box: 1