


## SHIMADZU Prominence HPLC System



Complete HPLC System with manual injection

Rating: Not Rated Yet

**Price**

Sales price \$18,250.00

Tax amount

[Ask a question about this product](#)

Manufacturer [Shimadzu](#)

### Description

CBM-20A	System Controller
DGU-20A5R	Degasser: 5 lines
LC-20AT	Solvent Delivery System: Low or High-pressure gradient modes
LC-20AD/AT	Low Pressure Gradient Unit: Combines up to four solvents
CTO-20AC	Column Heater: 10 °C below room temperature to 85 °C
SPD-20A	UV Detector: Single wavelength, dual wavelength and wavelength scanning

Rheodyne Model 7725i Manual Injector

Includes:

Solvent Tray

Software Version 5.57 SP1

Laptop Dell LATTITUDE 3440

Licenses and manuals

Solvent containers