


SHIMADZU Prominence HPLC System



Complete HPLC System with manual injection

Rating: Not Rated Yet

Price

Sales price \$18,250.00

Tax amount

[Ask a question about this product](#)

Manufacturer [Shimadzu](#)

Description

| | |
|------------|---|
| CBM-20A | System Controller |
| DGU-20A5R | Degasser: 5 lines |
| LC-20AT | Solvent Delivery System: Low or High-pressure gradient modes |
| LC-20AD/AT | Low Pressure Gradient Unit: Combines up to four solvents |
| CTO-20AC | Column Heater: 10 °C below room temperature to 85 °C |
| SPD-20A | UV Detector: Single wavelength, dual wavelength and wavelength scanning |

Rheodyne Model 7725i Manual Injector

Includes:

Solvent Tray

Software Version 5.57 SP1

Laptop Dell LATTITUDE 3440

Licenses and manuals

Solvent containers