

## SP Scientific Advantage Pro EL-85 Freeze Dryer / Lyophilizer with Leybold Trivac D4B Vacuum Pump



Full color touchscreen display

Rating: Not Rated Yet

**Price**

Sales price \$34,495.00

Tax amount

[Ask a question about this product](#)

Manufacturer [SP Scientific](#)

### Description

Flawless, Like New Condition

Uses Intellitronics™ controller

Full color touchscreen display

Three drying shelves

Condenser capacity: 6 Liters

Four thermal probes for direct measurement of samples

Minimum temperature -55 °C

Includes stainless steel wheeled cart

Includes Leybold Trivac D4B Vacuum Pump

Vacuum Pump Rate: 4 cubic meters per hour

Vacuum Pump Year of Manufacture: 2014

240 Volts