

## Spectronics Spectroline CM-10A Fluorescence Analysis Cabinet with EA-160 UV Light



UV Viewing Box

Rating: Not Rated Yet

**Price**

Sales price \$195.00

Tax amount

[Ask a question about this product](#)

Manufacturer [Spectronics Corporation](#)

### Description

UV Viewing Box

- Black fabric screen for easy hand insertion of objects;
- Floor can be easily removed for larger items;
- Front insertion panel dimensions: 10.75" W x 5.0" H
- Internal box depth is 8.50 inches;
- 120 Volts AC.