

Thermo Scientific Model 1387 1300 Series A2 6' Biological Safety Cabinet with Stand



Class II Type A2

Rating: Not Rated Yet

Price

Sales price \$7,245.00

Tax amount

[Ask a question about this product](#)

Manufacturer [Thermo Scientific](#)

Description

Thermo Scientific Model 1387 1300 Series A2 6' Biological Safety Cabinet with Stand
Like New in Every Way
Six Foot with stand

Series meets NSF/ANSI™ Standard 49, NSF/ANSI 49, UL, CE. Real-time airflow adjustments ensure inflow and downflow velocities remain steady.

Pressure sensor filter monitoring measures total airflow across the working area.

Proprietary SmartClean front window design simplifies cleaning for a safer working environment.

SmartFlow Indicator adjusts airflows as filter resistance changes.

SmartFlow visually demonstrates cabinet's ability to compensate to maintain personal protection.