

## US Scientific USDW-HL678S -86 C Freezer [New]



23.7 Cubic Feet

Rating: Not Rated Yet

**Price**

Sales price \$9,388.00

Tax amount

[Ask a question about this product](#)

Manufacturer [US Scientific](#)

### Description

US Scientific USDW-HL678S -86 C Freezer:

- Brand NEW: Has never been used
- Capacity: 23.7 Cu. Ft.
- Lockable outer door
- Uses single compressor
- Stainless steel interior
- Battery backup for controller
- Two insulated interior doors
- Three stainless steel shelves
- 7 inch LED touch screen controller
- USB port for retrieving stored data
- Password can be set at the control panel
- Includes 4 industrial casters & 4 leveling feet
- Operates on 115 VAC, 60 Hz main power (12.8 Amps)
- Includes a 1-year factory warranty from the manufacturer

