


BIOBASE BKM-Z80BIII B Series Table Top Autoclave, New



Capacity 2.8 Cubic Feet, One year warranty.

Rating: Not Rated Yet

Price

Sales price \$11,995.00

Tax amount

[Ask a question about this product](#)

Manufacturer [BIOBASE](#)

Description

BIOBASE BKM-Z80BIII B Series Table Top Autoclave, New:

- Brand NEW - never used!
- LCD display
- Internal water filter
- Water quality monitor
- Door interlock for safety
- Design temperature 144° C
- Temperature range 105° C > 136° C
- Vacuum system operates quietly
- Capacity 2.8 Cubic Feet (80 Liters)
- Build-in easy to clean open type water tank
- Medical-grade 304 Stainless Steel (SST) chamber
- Dual chamber overpressure protection for safety

- Operates on 208 VAC, 60 Hz, 3-phase main power
- Date of Manufacture: October 31, 2017
- Includes 1-year warranty directly from the manufacturer & backed directly by A.I.