

## CEM Smart 5 Microwave Moisture Analyzer



Uses microwave heat to measure moisture content

Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)

Manufacturer [CEM](#)

### Description

Uses microwave heat to measure moisture content

Use of microwave heat speeds up the measurement process

LCD display

Large internal working volume

Integrated printer for results

Able to control sample to specified temperature to prevent degradation

Maximum Power 300 Watts

Power output can be set lower for finer control of heating

Balance precision is 0.0001 grams

Integrated software can store up to 100 methods

120 Volts

