

## CEM Smart 5 Microwave Moisture Analyzer

Uses microwave heat to measure moisture content

Rating: Not Rated Yet

**Price**

Sales price \$14,485.00

Tax amount

[Ask a question about this product](#)

Manufacturer [CEM](#)

### Description

CEM Smart 5 Microwave Moisture Analyzer:

- Maximum Power 300 Watts
- Uses color LCD front panel screen
- Integrated printer for results
- Large internal working volume
- Balance precision is 0.0001 grams
- Integrated software can store up to 100 methods
- Uses microwave heat to measure moisture content
- Use of microwave heat speeds up the measurement process
- Able to control sample to specified temperature to prevent degradation
- Power output can be set lower for finer control of heating
- Operates on 120 VAC, 60 Hz, single phase main power
- Includes 90-day warranty directly from A.I.