

## Fisher Scientific 945415 Digital Vortex Mixer



Countdown timer for consistent results

Rating: Not Rated Yet

**Price**

Sales price \$95.00

Tax amount

[Ask a question about this product](#)

Manufacturer [Fisher Scientific](#)

### Description

Fisher Scientific Digital Vortex Mixer (945415):

- Countdown timer for consistent results
- Continuous duty or touch modes
- LED display for speed and time
- Speed Range: 500 RPM to 3,000 RPM
- Requires 115 VAC, 50/60 Hz main power
- Vortex stops when countdown timer reaches zero
- Includes 90-day warranty