

Fisher Scientific Isotemp 145D Dry Bath Block Heater



Digital Control

Rating: Not Rated Yet

Price

Sales price \$1,185.00

Tax amount

[Ask a question about this product](#)

Manufacturer [Fisher Scientific](#)

Description

Fisher Scientific Isotemp 145D Dry Bath Block Heater:

- Digital Control
- Bright, easy to read LED digital display screen
- Holds four (4) blocks
- Temperature Range: 5° C above ambient to a max of 130° C
- Steel body is powder coated and corrosion-resistant
- The aluminum heat blocks are highly conductive
- Operates on 115 VAC, 60 Hz, single-phase main power
- Includes 90-day warranty directly from A.I.