

HP 6890 Series GC System



Highly effective for detection of volatile compounds

Rating: Not Rated Yet

Price

Sales price \$37,000.00

Tax amount

[Ask a question about this product](#)

Manufacturer [HP](#)

Description

HP (Agilent) 6890 Series GC System / 5973 MS Turbo:

- Includes dual split / splitless injector, TCD, FID & Mass Spectrometer detectors
- NIST Library loaded on HP Compac computer running Chemstation (system given PM before shipping)
- CTC Analytics Combi Pal Headspace / Liquid Autosampler includes the following 5 items
- (...Allows user to use both liquid and headspace injections):
- (1) Liquid syringe kit,
- (2) Wash station,
- (3) Agitator/ incubator,
- (4) Tray set vial holders to be determined, and
- (5) Software integrating with Chemstation.

System also includes 90-day warranty directly from A.I.