

Mettler Model AE 240-S SNR L45095 Tabletop Balance



Mettler Model AE 240-S Tabletop Balance

Rating: Not Rated Yet

Price

Sales price \$1,500.00

Tax amount

[Ask a question about this product](#)

Manufacturer [Mettler](#)

Description

Mettler Model AE 240-S Tabletop Balance:

Description:

- 40g/200g
- Readability: 0.01mg/0.1mg
- Weighing Range: 0...41g/0...205g

- Tare Range: 0...41g/0...205g
- Reproducibility: 0.02mg/0.1mg
- Stabilizing Time: 8sec/5sec
- Includes 90-day warranty directly from A.I.