

## Mettler Toledo HB43 Moisture Analyzer



Halogen Moisture Balance

Rating: Not Rated Yet

**Price**

Sales price \$3,200.00

Tax amount

[Ask a question about this product](#)

Manufacturer [Mettler](#)

### Description

Mettler Toledo HB43 Moisture Analyzer.

- Adjustable drying temperature
- Includes black 6' AC power cord
- Uses Halogen bulb to heat sample
- Available weighing space 12 inch diameter
- Operates on 120 VAC, 60 Hz, single phase main power
- Includes 90-day warranty directly from A.I.