



Persee T8DS Double Beam UV-Vis Spectrophotometer [New]



Double Beam

Price with Rebate \$6,599.00

Discount

[Ask a question about this product](#)

Manufacturer [Persee](#)

Description

T8DS is a high performance double beam scanning spectrophotometer with a variable spectral

bandwidth.

The Czerny-Turner Monochromator with a holographic grating keeps stray light at a minimum and offers excellent optical resolution.

- The use of a photomultiplier tube as detector delivers exceptional sensitivity.
 - The true double-beam optical system coupled with an efficient and well proven electronic control system ensures high stability and low background noise.
 - High performance variable (0.1, 0.2, 0.5, 1.0, 2.0, 5.0nm) spectral bandwidth
 - Wavelength: 190-900nm
-