

## VWR SCUCBI-0404 Under Counter Refrigerator [Clearance]



Front ventilation design for placement between cabinetry

Price Not Set \$290.43

Discount

[Ask a question about this product](#)

Manufacturer [VWR](#)

Desktop Built-In Undercounter Laboratory Refrigerators feature a front ventilation design for built-in applications between cabinetry.

- Digital temperature display for convenient monitoring
- Remote alarm contacts
- Audible and visual alarm package for sample security

Standard features include forced draft circulation and a microprocessor digital temperature controller for superior temperature performance, UL listed.

#### Built-In Undercounter Refrigerator

- 120V, 5A
  - 130 L (4.6 cu. ft)
  - 84.8H x 60.3W x 61D cm (333/8 x 233/4 x 241/2")
  - Solid Door
  - 0 to 10°C (32 to 50°F)
  - 100 kg (220.5 lbs.)
  - Serial # SYM-448536-1403
  - Current listed price includes the clearance discount
  - Clearance items come with a 30-day warranty
-