

## Sartorius B120S Precision Balance



Capacity 121g, Precision 0.0001g

Rating: Not Rated Yet

**Price**

Sales price \$485.00

Tax amount

[Ask a question about this product](#)

Manufacturer [Sartorius](#)

### Description

Sartorius B120S Precision Balance:

- Capacity 121g, Precision 0.0001g
- Digital Display
- Internal Calibration
- Weigh Pan Diameter: 9 cm
- Includes 5-sided, glass draft shield
- Requires 115 VAC, 60 Hz, single phase main power
- Includes 90-day warranty directly from A.I.