

Thermo IEC Centra GP8F Undercounter Floor Centrifuge



Maximum Speed 5,300 RPM

Rating: Not Rated Yet

Price

Sales price \$1,185.00

Tax amount

[Ask a question about this product](#)

Manufacturer [Thermo](#)

Description

Thermo IEC Centra GP8F Undercounter Floor Centrifuge:

- Front panel digital control
- Able to store 35 programs
- Maximum Speed: 5,300 RPM
- External Dimensions: 23.0" W x 29.0" D x 35.0" H
- A wide variety of rotors & accessories are available
- Mounted on stand with storage drawer & lockable wheels
- Set speeds = 500 to 5300 RPM in 10 or 100 RPM increments
- Adjustable speed, time, temperature, acceleration & deceleration
- Operates on 120 VAC, 50/60 Hz, single-phase main power
- Can also operate on 220 VAC, 50/60 Hz, single-phase main power
- Includes 90-day warranty directly from A.I.