

Thermo Scientific 3566A Explosion-Proof Laboratory Refrigerator



Never used, 21 Cubic Feet

Rating: Not Rated Yet

Price

Sales price \$1,995.00

Tax amount

[Ask a question about this product](#)

Manufacturer [Thermo Scientific](#)

Description

Thermo Scientific 3566A Explosion-Proof Laboratory Refrigerator:

- Never used, missing kick plate and small exterior dents
- Lockable door
- Capacity is 21 cu. ft.
- Operates on 115 VAC, 60 Hz main power
- 4 fixed shelves, with integrated cooling
- Exterior Dimensions: 74.5" H x 32.0" W x 30.5" D
- Explosion proof for storage of highly flammable material
- Includes A.I.'s own 90-day warranty