

USA Scientific Model SD Micro Centrifuge



Speed: 6,000 RPM

Rating: Not Rated Yet

Price

Sales price \$125.50

Tax amount

[Ask a question about this product](#)

Manufacturer [USA Scientific](#)

Description

USA Scientific Model SD Micro Centrifuge:

- Fixed-angle rotor
- Speed: 6,000 RPM, 2000 x g
- 6 places \ positions (for 3 sizes of tubes)
- Holds 1.5 or 2.0 ml tubes (or use adapters for 0.5 ml tubes)
- Stops rotation with On\Off switch on side or by opening the lid
- Dimensions: 5.125" W x 5.875" D x 4.625" H.
- Weight: 2.27 lbs.
- Like new condition
- Operates on 115 VAC, 60 Hz main power