


## VWR 200 Mini Stirrer



Ideal for small lab spaces

Rating: Not Rated Yet

**Price**

Sales price \$87.50

Tax amount

[Ask a question about this product](#)

Manufacturer [VWR](#)

### Description

VWR 200 Mini Stirrer:

- Like new condition
- Durable cast aluminum surface
- Work surface: ~4 inches x 5 inches
- Includes built-in back-mounted ring stand clamp
- Analog speed control - ranges from 100 to 1200 RPM
- Operates on 115 VAC, 60 Hz main power
- Weight: 2.32 lbs.
- Dimensions: 3.375" H x 4.0" W x 6.25" D (includes built-in clamp on back).