

VWR 2350MP Water Jacketed Double Stack CO2 Incubator



Water Jacketed for temperature stability

Rating: Not Rated Yet

Price

Sales price \$3,645.00

Tax amount

[Ask a question about this product](#)

Manufacturer [VWR](#)

Description

VWR 2350MP Water Jacketed Double Stack CO₂ Incubator:

- Stainless Steel (SST) interior
- Water Jacketed for temperature stability
- Digital Control for temperature and CO₂ percentage
- Operates on 120 VAC, 60 Hz single-phase main power
- Includes 90-day warranty directly from A.I.