

## VWR Sheldon 2300 Water Jacketed CO2 Incubator



Water Jacketed for temperature stability

Rating: Not Rated Yet

**Price**

Sales price \$995.00

Tax amount

[Ask a question about this product](#)

Manufacturer [VWR](#)

### Description

VWR Sheldon 2300 Water Jacketed CO2 Incubator:

- Includes Digital Display, Digital Control
- Water Jacketed for temperature stability
- Exterior Dimensions: 25.0" W x 39.0" H x 21.5" D
- Operates on 120 VAC, 60 Hz, single phase main power